

## Product datasheet for **SP2019HRP**

### Chicken IgM Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgM Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	<b>ELISA</b> (1/10,000-1/100,000).
Reactivity:	Chicken
Host:	Goat
Immunogen:	Purified chicken IgM.
Formulation:	PBS containing 0.1% Proclin 300 as preservative and 0.2% BSA as stabilizer Label: HRP State: Liquid purified Ig fraction Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography.
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)