

## Product datasheet for **SP2019P**

### Chicken IgM Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgM Goat Polyclonal Antibody
Applications:	ELISA, ID
Recommended Dilution:	<b>ELISA:</b> (1/100-1/1000). <b>Immunodiffusion.</b>
Reactivity:	Chicken
Host:	Goat
Immunogen:	Purified chicken IgM.
Formulation:	PBS, pH 7.2 containing 0.09% Sodium Azide as preservative. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)