

Product datasheet for **SP2020HRP**

Chicken IgA Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgA Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA (1/10,000-1/100,000).
Reactivity:	Chicken
Host:	Goat
Immunogen:	Purified chicken IgA.
Formulation:	PBS containing 0.1% Proclin 300 as preservative and 0.2% BSA as stabilizer. Label: HRP State: Liquid purified IgG fraction. Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography.
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)