

Product datasheet for **SP2021F**

Chicken IgG / Chicken IgY (Fc specific) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgG / Chicken IgY (Fc specific) Goat Polyclonal Antibody
Applications:	FC, IHC
Recommended Dilution:	Flow cytometry.
Reactivity:	Chicken
Host:	Goat
Immunogen:	Purified chicken IgG - Fc fragment
Formulation:	PBS containing 0.09% Sodium Azide and 0.2% Bovine Serum Albumin Label: FITC State: Liquid purified IgG Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light.



[View online »](#)