

## Product datasheet for **SP2021HRP**

### Chicken IgG / Chicken IgY (Fc specific) Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgG / Chicken IgY (Fc specific) Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 1/10,000-1/100,000.
Reactivity:	Chicken
Host:	Goat
Immunogen:	Purified Chicken IgG (Fc) fragment.
Formulation:	PBS containing 0.1% ProClin™ 300 as preservative and 0.2% BSA as stabilizer. Label: HRP State: Liquid purified IgG fraction. Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C. <b>DO NOT FREEZE!</b> Should this product contain a precipitate we recommend microcentrifugation before use.
Note:	Antiserum Preparation: Antisera to Chicken IgG were raised by repeated immunisation of goat with highly purified antigen.



[View online »](#)