

Product datasheet for **SP2021P**

Chicken IgG / Chicken IgY (Fc specific) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgG / Chicken IgY (Fc specific) Goat Polyclonal Antibody
Applications:	ELISA, ID
Recommended Dilution:	ELISA (1/100-1/10000). Immunodiffusion.
Reactivity:	Chicken
Host:	Goat
Immunogen:	Purified chicken IgG - Fc fragment
Formulation:	PBS, pH 7.2 containing 0.09% Sodium Azide as preservative. State: Purified State: Liquid purified IgG fraction.
Concentration:	lot specific
Purification:	Affinity Chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)