

Product datasheet for **SP2041B**

Rabbit IgG (H+L chain) Sheep Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit IgG (H+L chain) Sheep Polyclonal Antibody
Applications:	ELISA, IHC
Recommended Dilution:	ELISA: 1/25600 - 1/51200. Western Blot: 1/20000. Immunohistochemistry on frozen sections: 1/400 - 1/800. Immunohistochemistry on paraffin sections: 1/100 - 1/200.
Reactivity:	Rabbit
Host:	Sheep
Immunogen:	Rabbit IgG from pooled sera
Formulation:	PBS, pH 7.2 containing 0.09% Sodium Azide and 1% Bovine Serum Albumin Label: Biotin State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)