

## Product datasheet for SP2083F

## Guinea pig IgG (H+L chain) Goat Polyclonal Antibody

## **Product data:**

Product Type:	Secondary Antibodies
Product Name:	Guinea pig IgG (H+L chain) Goat Polyclonal Antibody
Applications:	FC, IHC
Recommended Dilution:	Flow Cytometry. Immunohistochemistry on frozen sections: 1:200.
Reactivity:	Guinea Pig
Host:	Goat
Formulation:	PBS, pH 7.2 containing 0.09% Sodium Azide and 0.2% Bovine Serum Albumin Label: FITC State: Liquid purified IgG Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light.

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US