

## Product datasheet for **SP2083HRP**

### Guinea pig IgG (H+L chain) Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Guinea pig IgG (H+L chain) Goat Polyclonal Antibody
Applications:	ELISA, IHC
Recommended Dilution:	<b>ELISA</b> (1/1000-1/10,000). <b>Immunohistochemistry on Frozen Sections.</b> <b>Immunohistochemistry on Paraffin Sections.</b>
Reactivity:	Guinea Pig
Host:	Goat
Immunogen:	Highly purified Guinea pig IgG (H&L)
Isotype:	IgG
Formulation:	PBS Label: HRP State: Purified State: Liquid purified IgG fraction Stabilizer: 0.2% BSA Preservative: 0.1% Proclin 300™ Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C. <b>DO NOT FREEZE!</b>



[View online »](#)