

Product datasheet for SP2086F

Porcine IgG (H+L chain) Goat Polyclonal Antibody

Product data:

Product Type: Secondary Antibodies

Product Name: Porcine IgG (H+L chain) Goat Polyclonal Antibody

Applications: FC, IHC

Recommended Dilution: Flow Cytometry.

Immunohistochemistry on frozen sections: 1:200 - 1:2000.

Reactivity: Porcine Host: Goat

Formulation: PBS, pH 7.2 containing 0.09% Sodium Azide and 0.2% BSA

Label: FITC

State: Liquid purified IgG

Label: Fluorescein Isothiocyanate Isomer 1

Concentration: lot specific

Purification: Affinity chromatography

Conjugation: FITC

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

This product is photosensitive and should be protected from light.

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com