

Product datasheet for **SP2086HRP**

Porcine IgG (H+L chain) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Porcine IgG (H+L chain) Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1:10,000 - 1:100,000.
Reactivity:	Porcine
Host:	Goat
Formulation:	PBS containing 0.09% Proclin (TM)300 as preservative and 0.2% Bovine Serum Albumin as stabilizer Label: HRP State: Liquid purified IgG Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)