

Product datasheet for **SP2124HRP**

Canine IgE Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Canine IgE Goat Polyclonal Antibody
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA: 1/10,000-1/100,000. Western Blot: 1/1000-1/10,000. Immunohistochemistry on Frozen Sections: 1/200-1/500.
Reactivity:	Canine
Host:	Goat
Immunogen:	Highly purified Canine IgE
Formulation:	PBS Label: HRP State: Liquid purified IgG fraction Stabilizer: 0.2% BSA Preservative: 0.1% Proclin™ 300 Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)