

## Product datasheet for **SP2124P**

### Canine IgE Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Canine IgE Goat Polyclonal Antibody
Applications:	ELISA, FC, WB
Recommended Dilution:	<b>ELISA:</b> 1/100-1/10,000. <b>Western Blot:</b> 1/1000-1/10,000. <b>Flow Cytometry:</b> 1/200-1/2,000. This antibody has been reported for use in <b>Immunohistology of Frozen Canine Sections</b> .
Reactivity:	Canine
Host:	Goat
Immunogen:	Highly purified Canine IgE.
Formulation:	PBS, pH 7.2 containing 0.09% Sodium Azide as preservative State: Purified State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Affinity Chromatography.
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)