

## Product datasheet for **SP2125HRP**

### Canine IgG1 Goat Polyclonal Antibody

#### Product data:

|                       |  |
|-----------------------|--|
| Product Type:         | Secondary Antibodies   |
| Product Name:         | Canine IgG1 Goat Polyclonal Antibody   |
| Applications:         | ELISA, IHC, WB   |
| Recommended Dilution: | ELISA: 1:10,000 - 1:100,000.<br>Western Blot: 1:1000 - 1:10,000.<br>Immunohistochemistry on frozen sections: 1:200 - 1:500.                              |
| Reactivity:           | Canine   |
| Host:                 | Goat   |
| Formulation:          | PBS containing 0.1% Proclin <sup>TM</sup> 300 and 0.2% Bovine Serum Albumin<br>Label: HRP<br>State: Liquid purified IgG<br>Label: Horseradish Peroxidase |
| Concentration:        | lot specific   |
| Purification:         | Affinity chromatography  |
| Conjugation:          | HRP  |
| Storage:              | Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.<br>Avoid repeated freezing and thawing.                                     |



[View online »](#)