

## Product datasheet for **SP2125P**

### Canine IgG1 Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Canine IgG1 Goat Polyclonal Antibody
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA: 1:100 - 1:10,000. Western Blot: 1:1000 - 1:10,000. Immunohistochemistry on frozen sections: 1:200 - 1:2000.
Reactivity:	Canine
Host:	Goat
Formulation:	PBS containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)