

Product datasheet for **SP2130F**

Sheep IgM Rabbit Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Sheep IgM Rabbit Polyclonal Antibody
Applications:	FC
Recommended Dilution:	Flow Cytometry: 1:200 - 1:2000.
Reactivity:	Sheep
Host:	Rabbit
Formulation:	PBS, pH 7.2 containing 0.09% Sodium Azide and 0.2% Bovine Serum Albumin Label: FITC State: Liquid purified IgG Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light.



[View online »](#)