

## Product datasheet for **SP2130HRP**

### Sheep IgM Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Sheep IgM Rabbit Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1/10,000-1/100,000.
Reactivity:	Sheep
Host:	Rabbit
Formulation:	PBS containing 0.1% ProClin300 as preservative and 0.2% BSA as stabilizer. Label: HRP State: Liquid purified IgG fraction. Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography.
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)