

Product datasheet for SP2130P

Sheep IgM Rabbit Polyclonal Antibody

Product data:

Product Type: Secondary Antibodies

Product Name: Sheep IgM Rabbit Polyclonal Antibody

Applications: ELISA

Recommended Dilution: ELISA: 1:100 - 1:10,000.

Reactivity: Sheep Host: Rabbit

Formulation: PBS, pH 7.2 containing 0.09% Sodium Azide

State: Purified

State: Liquid purified IgG

Concentration: lot specific

Purification: Affinity chromatography

Conjugation: Unconjugated

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com