

Product datasheet for **SP2131HRP**

Canine IgG2 Sheep Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Canine IgG2 Sheep Polyclonal Antibody
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA: 1:10,000 - 1:100,000. Western Blot: 1:1000 - 1:10,000. Immunohistochemistry on frozen sections: 1:200 - 1:500.
Reactivity:	Canine
Host:	Sheep
Formulation:	PBS containing 0.09% Proclin™ ₃₀₀ and 0.2% BSA State: Liquid purified IgG conjugated to HRP Label: Horseradish Peroxidase (HRP)
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE. Avoid repeated freezing and thawing.



[View online »](#)