

Product datasheet for **SP2185AP**

Mouse lambda light chain Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse lambda light chain Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA (1/1,000-1/2,000).
Reactivity:	Mouse
Host:	Goat
Immunogen:	Mouse paraproteins with lambda 1 and lambda 2 light chains
Formulation:	50mM TRIS, 1mM MgCl ₂ , 0.09% Sodium Azide (NaN ₃), 50% Glycerol Label: AP State: Liquid Label: Alkaline Phosphatase
Purification:	Purified Ig was prepared from whole serum by affinity chromatography on pooled mouse Igs with lambda chains covalently linked to agarose.
Conjugation:	AP
Storage:	Store at 2-8°C or at -20°C if preferred. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use. Shelf life: One year from despatch.
Note:	Store at 2-8°C or at -20°C if preferred. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use. Shelf life: One year from despatch.



[View online »](#)