

Product datasheet for **SP2185B**

Mouse lambda light chain Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse lambda light chain Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1/1000 - 1/2000.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Mouse paraproteins with lambda 1 and lambda 2 light chains
Formulation:	PBS containing 0.09% Sodium Azide Label: Biotin State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)