

Product datasheet for **SP2185F**

Mouse lambda light chain Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse lambda light chain Goat Polyclonal Antibody
Applications:	FC, IF
Recommended Dilution:	Flow cytometry. Immunofluorescence.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Mouse paraproteins with lambda 1 and lambda 2 light chains
Formulation:	PBS containing 0.09% sodium azide as preservative Label: FITC State: Liquid IgG fraction Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity chromatography on pooled mouse Igs with lambda chains covalently linked to agarose
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light.



[View online »](#)