

Product datasheet for SP2185HRP

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn

OriGene Technologies, Inc.

Mouse lambda light chain Goat Polyclonal Antibody

Product data:

Product Type: Secondary Antibodies

Product Name: Mouse lambda light chain Goat Polyclonal Antibody

Applications: ELISA

Recommended Dilution: ELISA (1/4,000-1/8,000).

Reactivity: Mouse **Host:** Goat

Immunogen: Mouse paraproteins with lambda 1 and lambda 2 light chains

Formulation: Phosphate buffered saline, 50% Glycerol

Label: HRP

State: Liquid Ig fraction

Label: Horseradish Peroxidase

Purification: Affinity chromatography on pooled mouse Igs with lambda chains covalently linked to

agarose

Conjugation: HRP

Storage: Store the antibody at 2-8 °C for up to one week and at -20°C for longer. This product should

be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.