

## Product datasheet for **SP2185P**

### Mouse lambda light chain Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse lambda light chain Goat Polyclonal Antibody
Applications:	ELISA, IF, WB
Recommended Dilution:	ELISA: 1/4000 - 1/8000. Western Blot. Immunofluorescence. Immunoblotting.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Mouse paraproteins with lambda 1 and lambda 2 light chains
Formulation:	Borate buffered saline State: Aff - Purified State: Liquid Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.



[View online »](#)