

## Product datasheet for **SR301340**

### EEF1G Human siRNA Oligo Duplex (Locus ID 1937)

#### Product data:

|                     |                                                                                                                                                                                                                                                                                                                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:       | siRNA Oligo Duplexes                                                                                                                                                                                                                                                                                                                                                                                |
| Purity:             | HPLC purified                                                                                                                                                                                                                                                                                                                                                                                       |
| Quality Control:    | Tested by ESI-MS                                                                                                                                                                                                                                                                                                                                                                                    |
| Sequences:          | Available with shipment                                                                                                                                                                                                                                                                                                                                                                             |
| Stability:          | One year from date of shipment when stored at -20°C.                                                                                                                                                                                                                                                                                                                                                |
| # of transfections: | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).                                                                                                                                                                                                                                                                                              |
| Note:               | Single siRNA duplex (10nmol) can be ordered.                                                                                                                                                                                                                                                                                                                                                        |
| RefSeq:             | <a href="#">NM_001404</a>                                                                                                                                                                                                                                                                                                                                                                           |
| UniProt ID:         | <a href="#">P26641</a>                                                                                                                                                                                                                                                                                                                                                                              |
| Synonyms:           | EF1G; GIG35                                                                                                                                                                                                                                                                                                                                                                                         |
| Components:         | EEF1G (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 1937)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml                                                                                                                                                |
| Summary:            | This gene encodes a subunit of the elongation factor-1 complex, which is responsible for the enzymatic delivery of aminoacyl tRNAs to the ribosome. This subunit contains an N-terminal glutathione transferase domain, which may be involved in regulating the assembly of multisubunit complexes containing this elongation factor and aminoacyl-tRNA synthetases. [provided by RefSeq, Jul 2008] |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).