

## Product datasheet for **SR302197**

### **HOXB8 Human siRNA Oligo Duplex (Locus ID 3218)**

#### **Product data:**

|                            |   |
|----------------------------|---|
| <b>Product Type:</b>       | siRNA Oligo Duplexes  |
| <b>Purity:</b>             | HPLC purified   |
| <b>Quality Control:</b>    | Tested by ESI-MS  |
| <b>Sequences:</b>          | Available with shipment   |
| <b>Stability:</b>          | One year from date of shipment when stored at -20°C.  |
| <b># of transfections:</b> | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).  |
| <b>Note:</b>               | Single siRNA duplex (10nmol) can be ordered.  |
| <b>RefSeq:</b>             | <a href="#">NM_024016</a>   |
| <b>UniProt ID:</b>         | <a href="#">P17481</a>  |
| <b>Synonyms:</b>           | Hox-2.4; HOX2; HOX2D  |
| <b>Components:</b>         | HOXB8 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 3218)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml  |
| <b>Summary:</b>            | This gene is a member of the Antp homeobox family and encodes a nuclear protein with a homeobox DNA-binding domain. It is included in a cluster of homeobox B genes located on chromosome 17. The encoded protein functions as a sequence-specific transcription factor that is involved in development. Increased expression of this gene is associated with colorectal cancer. Mice that have had the murine ortholog of this gene knocked out exhibit an excessive pathologic grooming behavior. This behavior is similar to the behavior of humans suffering from the obsessive-compulsive spectrum disorder trichotillomania. [provided by RefSeq, Jul 2008] |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).