

## Product datasheet for **SR302572**

### **KIR2DS4 Human siRNA Oligo Duplex (Locus ID 3809)**

#### **Product data:**

|                            |   |
|----------------------------|---|
| <b>Product Type:</b>       | siRNA Oligo Duplexes  |
| <b>Purity:</b>             | HPLC purified   |
| <b>Quality Control:</b>    | Tested by ESI-MS  |
| <b>Sequences:</b>          | Available with shipment   |
| <b>Stability:</b>          | One year from date of shipment when stored at -20°C.  |
| <b># of transfections:</b> | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).  |
| <b>Note:</b>               | Single siRNA duplex (10nmol) can be ordered.  |
| <b>RefSeq:</b>             | <a href="#">NM_001281971</a> , <a href="#">NM_001281972</a> , <a href="#">NM_012314</a> , <a href="#">NM_178228</a>   |
| <b>UniProt ID:</b>         | <a href="#">P43632</a>  |
| <b>Synonyms:</b>           | CD158I; KIR-2DS4; KIR1D; KIR412; KKA3; NKAT-8; NKAT8  |
| <b>Components:</b>         | KIR2DS4 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 3809)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml  |
| <b>Summary:</b>            | Killer cell immunoglobulin-like receptors (KIRs) are transmembrane glycoproteins expressed by natural killer cells and subsets of T cells. The KIR genes are polymorphic and highly homologous and they are found in a cluster on chromosome 19q13.4 within the 1 Mb leukocyte receptor complex (LRC). The gene content of the KIR gene cluster varies among haplotypes, although several "framework" genes are found in all haplotypes (KIR3DL3, KIR3DP1, KIR3DL4, KIR3DL2). The KIR proteins are classified by the number of extracellular immunoglobulin domains (2D or 3D) and by whether they have a long (L) or short (S) cytoplasmic domain. KIR proteins with the long cytoplasmic domain transduce inhibitory signals upon ligand binding via an immune tyrosine-based inhibitory motif (ITIM), while KIR proteins with the short cytoplasmic domain lack the ITIM motif and instead associate with the TYRO protein tyrosine kinase binding protein to transduce activating signals. The ligands for several KIR proteins are subsets of HLA class I molecules; thus, KIR proteins are thought to play an important role in regulation of the immune response. [provided by RefSeq, Jul 2008] |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).