

## Product datasheet for **SR302593**

### SRP1 (KPNA1) Human siRNA Oligo Duplex (Locus ID 3836)

#### Product data:

|                     |  |
|---------------------|--|
| Product Type:       | siRNA Oligo Duplexes   |
| Purity:             | HPLC purified  |
| Quality Control:    | Tested by ESI-MS   |
| Sequences:          | Available with shipment  |
| Stability:          | One year from date of shipment when stored at -20°C.   |
| # of transfections: | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).   |
| Note:               | Single siRNA duplex (10nmol) can be ordered.   |
| RefSeq:             | <a href="#">NM_002264</a> , <a href="#">NR_026698</a>  |
| UniProt ID:         | <a href="#">P52294</a>   |
| Synonyms:           | IPOA5; NPI-1; RCH2; SRP1   |
| Components:         | KPNA1 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 3836)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml   |
| Summary:            | The transport of molecules between the nucleus and the cytoplasm in eukaryotic cells is mediated by the nuclear pore complex (NPC), which consists of 60-100 proteins. Small molecules (up to 70 kD) can pass through the nuclear pore by nonselective diffusion while larger molecules are transported by an active process. The protein encoded by this gene belongs to the importin alpha family, and is involved in nuclear protein import. This protein interacts with the recombination activating gene 1 (RAG1) protein and is a putative substrate of the RAG1 ubiquitin ligase. Alternative splicing results in multiple transcript variants.<br>[provided by RefSeq, Nov 2012] |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).