

Product datasheet for **SR305212**

zinc finger protein 138 (ZNF138) Human siRNA Oligo Duplex (Locus ID 7697)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_001160183 , NM_001271637 , NM_001271638 , NM_001271639 , NM_001271640 , NM_001271649 , NM_006524 , NR_027660 , NR_027661 , NR_073389 , NR_073391 , NR_160269 , NM_001367573 , NM_001367575 , NM_001367577 , NR_160268 , NR_160270 , NR_160271 , NM_001367572 , NM_001367574 , NM_001367576
UniProt ID:	P52744
Synonyms:	pHZ-32
Components:	ZNF138 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 7697) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	May be involved in transcriptional regulation as a repressor.[UniProtKB/Swiss-Prot Function]
Performance Guaranteed:	OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).



[View online »](#)