

Product datasheet for **SR305247**

ZNF208 Human siRNA Oligo Duplex (Locus ID 7757)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_001329971 , NM_001329972 , NM_001329973 , NM_001329974 , NM_007153 , NR_138252
UniProt ID:	Q43345
Synonyms:	PMIDP; ZNF95
Components:	ZNF208 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 7757) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Zinc finger proteins (ZNFs), such as ZNF208, bind DNA and, through this binding, regulate gene transcription. Most ZNFs contain conserved C2H2 motifs and are classified as Kruppel-type zinc fingers. A conserved protein motif, termed the Kruppel-associated box (KRAB) domain, mediates protein-protein interactions (Eichler et al., 1998 [PubMed 9724325]). See ZNF91 (MIM 603971) for further information on ZNFs.[supplied by OMIM, Aug 2009]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).