

## Product datasheet for **SR306839**

### ARFRP1 Human siRNA Oligo Duplex (Locus ID 10139)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_001134758</a> , <a href="#">NM_001267544</a> , <a href="#">NM_001267545</a> , <a href="#">NM_001267546</a> , <a href="#">NM_001267547</a> , <a href="#">NM_001267548</a> , <a href="#">NM_001267549</a> , <a href="#">NM_003224</a> , <a href="#">NR_051954</a> , <a href="#">NR_051955</a> , <a href="#">NR_051956</a> , <a href="#">NR_051957</a> , <a href="#">NR_051958</a>
UniProt ID:	<a href="#">Q13795</a>
Synonyms:	ARL18; ARP; Arp1
Components:	ARFRP1 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 10139) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNAse free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	The protein encoded by this gene is a membrane-associated GTP-ase which localizes to the plasma membrane and is related to the ADP-ribosylation factor (ARF) and ARF-like (ARL) proteins. This gene plays a role in membrane trafficking between the trans-Golgi network and endosomes. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, May 2012]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).