

Product datasheet for **SR306866**

SERF2 Human siRNA Oligo Duplex (Locus ID 10169)

Product data:

| | |
|---------------------|--|
| Product Type: | siRNA Oligo Duplexes |
| Purity: | HPLC purified |
| Quality Control: | Tested by ESI-MS |
| Sequences: | Available with shipment |
| Stability: | One year from date of shipment when stored at -20°C. |
| # of transfections: | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM). |
| Note: | Single siRNA duplex (10nmol) can be ordered. |
| RefSeq: | NM_001018108 , NM_001199875 , NM_001199876 , NM_001199877 , NM_001199878 , NM_005770 , NR_037672 |
| UniProt ID: | P84101 |
| Synonyms: | 4F5REL; FAM2C; H4F5REL; HsT17089 |
| Components: | SERF2 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 10169) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml |
| Summary: | Positive regulator of amyloid protein aggregation and proteotoxicity (PubMed:20723760). Induces conformational changes in amyloid proteins, such as HTT, driving them into compact formations preceding the formation of aggregates (PubMed:20723760).[UniProtKB/Swiss-Prot Function] |



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).