

## Product datasheet for **SR307971**

### NMNAT2 Human siRNA Oligo Duplex (Locus ID 23057)

#### Product data:

|                     |  |
|---------------------|--|
| Product Type:       | siRNA Oligo Duplexes   |
| Purity:             | HPLC purified  |
| Quality Control:    | Tested by ESI-MS   |
| Sequences:          | Available with shipment  |
| Stability:          | One year from date of shipment when stored at -20°C.   |
| # of transfections: | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).   |
| Note:               | Single siRNA duplex (10nmol) can be ordered.   |
| RefSeq:             | <a href="#">NM_015039</a> , <a href="#">NM_170706</a>  |
| UniProt ID:         | <a href="#">Q9BZQ4</a>   |
| Synonyms:           | C1orf15; PNAT2   |
| Components:         | NMNAT2 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 23057)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml   |
| Summary:            | This gene product belongs to the nicotinamide mononucleotide adenylyltransferase (NMNAT) enzyme family, members of which catalyze an essential step in NAD (NADP) biosynthetic pathway. Unlike the other human family member, which is localized to the nucleus, and is ubiquitously expressed; this enzyme is cytoplasmic, and is predominantly expressed in the brain. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008] |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).