

Product datasheet for **SR308080**

FBXL7 Human siRNA Oligo Duplex (Locus ID 23194)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_001278317 , NM_012304
UniProt ID:	Q9UJT9
Synonyms:	FBL6; FBL7
Components:	FBXL7 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 23194) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene encodes a member of the F-box protein family which is characterized by a 42-48 amino acid motif, the F-box, which binds to the S-phase kinase-associated protein 1 (Skp1) protein. The F-box proteins constitute one of the four subunits of E3 ubiquitin protein ligases called SCFs (SKP1-Cul1-F-box), which play a role in phosphorylation-dependent ubiquitination of proteins. The F-box proteins are divided into 3 subfamilies based on the other domain in the protein: F-box proteins that also have a WD-40 domain (Fbws subfamily), F-box proteins that also have leucine-rich repeats (Fbls subfamily) and F-box proteins that contain other motifs or lack known protein-interaction domains (Fbxs subfamily). The protein encoded by this gene belongs to the Fbls subfamily. Alternative splicing results in multiple transcript variants of this gene. [provided by RefSeq, May 2013]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).