

## Product datasheet for **SR308143**

### LPHN2 (ADGRL2) Human siRNA Oligo Duplex (Locus ID 23266)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_001297704</a> , <a href="#">NM_001297705</a> , <a href="#">NM_001297706</a> , <a href="#">NM_001330645</a> , <a href="#">NM_012302</a> , <a href="#">NM_001350698</a> , <a href="#">NM_001350699</a> , <a href="#">NM_001366002</a> , <a href="#">NM_001366004</a> , <a href="#">NM_001366006</a> , <a href="#">NM_001366008</a> , <a href="#">NM_001366003</a> , <a href="#">NM_001366005</a> , <a href="#">NM_001366007</a> , <a href="#">NM_001366009</a>
UniProt ID:	<a href="#">O95490</a>
Synonyms:	CIRL2; CL2; LEC1; LPHH1; LPHN2
Components:	ADGRL2 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 23266) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene encodes a member of the latrophilin subfamily of G-protein coupled receptors. The encoded protein participates in the regulation of exocytosis. The proprotein is thought to be further cleaved within a cysteine-rich G-protein-coupled receptor proteolysis site into two chains that are non-covalently bound at the cell membrane. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).