

## Product datasheet for **SR308420**

### **DNAAF11 Human siRNA Oligo Duplex (Locus ID 23639)**

#### **Product data:**

<b>Product Type:</b>	siRNA Oligo Duplexes
<b>Purity:</b>	HPLC purified
<b>Quality Control:</b>	Tested by ESI-MS
<b>Sequences:</b>	Available with shipment
<b>Stability:</b>	One year from date of shipment when stored at -20°C.
<b># of transfections:</b>	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
<b>Note:</b>	Single siRNA duplex (10nmol) can be ordered.
<b>RefSeq:</b>	<u><a href="#">NM_001321961</a></u> , <u><a href="#">NM_001321962</a></u> , <u><a href="#">NM_001321963</a></u> , <u><a href="#">NM_001321964</a></u> , <u><a href="#">NM_001321965</a></u> , <u><a href="#">NM_001321966</a></u> , <u><a href="#">NM_012472</a></u> , <u><a href="#">NR_073525</a></u> , <u><a href="#">NR_135905</a></u> , <u><a href="#">NR_135906</a></u> , <u><a href="#">NR_135907</a></u> , <u><a href="#">NR_135908</a></u> , <u><a href="#">NR_135909</a></u> , <u><a href="#">NR_135910</a></u> , <u><a href="#">NR_135911</a></u> , <u><a href="#">NR_135912</a></u> , <u><a href="#">NR_135913</a></u>
<b>UniProt ID:</b>	<u><a href="#">Q86X45</a></u>
<b>Synonyms:</b>	CILD19; LRRC6; LRTP; tiIB; TSLRP
<b>Components:</b>	LRRC6 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 23639) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
<b>Summary:</b>	The protein encoded by this gene contains several leucine-rich repeat domains and appears to be involved in the motility of cilia. Defects in this gene are a cause of primary ciliary dyskinesia-19 (CILD19). Alternative splicing of this gene results in multiple transcript variants. Related pseudogenes have been identified on chromosomes 4, 11 and 22. [provided by RefSeq, Apr 2016]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).