

## Product datasheet for **SR308743**

### **NAT9 Human siRNA Oligo Duplex (Locus ID 26151)**

#### **Product data:**

<b>Product Type:</b>	siRNA Oligo Duplexes
<b>Purity:</b>	HPLC purified
<b>Quality Control:</b>	Tested by ESI-MS
<b>Sequences:</b>	Available with shipment
<b>Stability:</b>	One year from date of shipment when stored at -20°C.
<b># of transfections:</b>	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
<b>Note:</b>	Single siRNA duplex (10nmol) can be ordered.
<b>RefSeq:</b>	<u><a href="#">NM_001305077</a></u> , <u><a href="#">NM_001305078</a></u> , <u><a href="#">NM_001305079</a></u> , <u><a href="#">NM_001305080</a></u> , <u><a href="#">NM_001305081</a></u> , <u><a href="#">NM_001305082</a></u> , <u><a href="#">NM_001305083</a></u> , <u><a href="#">NM_001305084</a></u> , <u><a href="#">NM_001305085</a></u> , <u><a href="#">NM_001305086</a></u> , <u><a href="#">NM_001305087</a></u> , <u><a href="#">NM_001305088</a></u> , <u><a href="#">NM_015654</a></u> , <u><a href="#">NR_130953</a></u> , <u><a href="#">NR_130954</a></u>
<b>UniProt ID:</b>	<u><a href="#">Q9BTE0</a></u>
<b>Synonyms:</b>	EBSP; hNATL
<b>Components:</b>	NAT9 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 26151) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNAse free siRNA Duplex Resuspension Buffer - 2 ml
<b>Performance Guaranteed:</b>	OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).



[View online »](#)