

## Product datasheet for **SR308876**

### **BSCL2 Human siRNA Oligo Duplex (Locus ID 26580)**

#### **Product data:**

|                            |  |
|----------------------------|--|
| <b>Product Type:</b>       | siRNA Oligo Duplexes   |
| <b>Purity:</b>             | HPLC purified  |
| <b>Quality Control:</b>    | Tested by ESI-MS   |
| <b>Sequences:</b>          | Available with shipment  |
| <b>Stability:</b>          | One year from date of shipment when stored at -20°C.   |
| <b># of transfections:</b> | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).   |
| <b>Note:</b>               | Single siRNA duplex (10nmol) can be ordered.   |
| <b>RefSeq:</b>             | <a href="#">NM_001122955</a> , <a href="#">NM_001130702</a> , <a href="#">NM_032667</a> , <a href="#">NR_037948</a> , <a href="#">NR_037949</a>  |
| <b>UniProt ID:</b>         | <a href="#">Q96G97</a>   |
| <b>Synonyms:</b>           | GNG3LG; HMN5; HMN5C; PELD; SPG17   |
| <b>Components:</b>         | BSCL2 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 26580)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml  |
| <b>Summary:</b>            | This gene encodes the multi-pass transmembrane protein protein seipin. This protein localizes to the endoplasmic reticulum and may be important for lipid droplet morphology. Mutations in this gene have been associated with congenital generalized lipodystrophy type 2 or Berardinelli-Seip syndrome, a rare autosomal recessive disease characterized by a near absence of adipose tissue and severe insulin resistance. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. Naturally occurring read-through transcription occurs between this locus and the neighboring locus HNRNPUL2 (heterogeneous nuclear ribonucleoprotein U-like 2).[provided by RefSeq, Mar 2011] |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).