

Product datasheet for **SR310403**

ULK4 Human siRNA Oligo Duplex (Locus ID 54986)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_001322500 , NM_001322501 , NM_017886 , NR_136342
UniProt ID:	Q96C45
Synonyms:	FAM7C1; REC01035
Components:	ULK4 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 54986) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene encodes a member of the unc-51-like serine/threonine kinase (STK) family. Members of this protein family play a role in neuronal growth and endocytosis. The encoded protein is likely involved in neurite branching, neurite elongation and neuronal migration. Genome-wide association studies (GWAS) indicate an association of variations in this gene with blood pressure and hypertension. Sequence variations in this gene may also be associated with psychiatric disorders, including schizophrenia and bipolar disorder. Pseudogenes associated with this gene have been identified and are located on chromosome 15. [provided by RefSeq, Jul 2016]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).