

Product datasheet for **SR310576**

ZNF532 Human siRNA Oligo Duplex (Locus ID 55205)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_001318726 , NM_001318727 , NM_001318728 , NM_018181 , NM_001353525 , NM_001353526 , NM_001353527 , NM_001353528 , NM_001353529 , NM_001353530 , NM_001353531 , NM_001353532 , NM_001353533 , NM_001353534 , NM_001353535 , NM_001353536 , NM_001353537 , NM_001353538 , NR_148459 , NR_148460
UniProt ID:	Q9HCE3
Components:	ZNF532 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 55205) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Performance Guaranteed:	OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).



[View online »](#)