

## Product datasheet for **SR311130**

### Leiomodin 3 (LMOD3) Human siRNA Oligo Duplex (Locus ID 56203)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<a href="#">NM_001304418</a> , <a href="#">NM_198271</a>
UniProt ID:	<a href="#">Q0VAK6</a>
Synonyms:	NEM10
Components:	LMOD3 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 56203) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	<p>The protein encoded by this gene is a member of the leiomodin family of proteins. This protein contains three actin-binding domains, a tropomyosin domain, a leucine-rich repeat domain, and a Wiskott-Aldrich syndrome protein homology 2 domain (WH2). Localization of this protein to the pointed ends of thin filaments has been observed, and there is evidence that this protein acts as a catalyst of actin nucleation, and is important to the organization of sarcomeric thin filaments in skeletal muscles. Mutations in this gene have been associated as one cause of Nemaline myopathy, as other genes have also been linked to this disorder. Nemaline myopathy is a disorder characterized by nonprogressive generalized muscle weakness and protein inclusions (nemaline bodies) in skeletal myofibers. Patients with mutations in this gene often present with a severe congenital form of the disorder. [provided by RefSeq, Jan 2015]</p>



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).