

Product datasheet for **SR311788**

HIVEP3 Human siRNA Oligo Duplex (Locus ID 59269)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_001127714 , NM_024503 , NR_038260 , NR_038261
UniProt ID:	Q5T1R4
Synonyms:	KBP-1; KBP1; KRC; Schnurri-3; SHN3; ZAS3; ZNF40C
Components:	HIVEP3 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 59269) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	This gene encodes a member of the human immunodeficiency virus type 1 enhancer-binding protein family. Members of this protein family contain multiple zinc finger and acid-rich (ZAS) domains and serine-threonine rich regions. This protein acts as a transcription factor and is able to regulate nuclear factor kappaB-mediated transcription by binding the kappaB motif in target genes. This protein also binds the recombination signal sequence that flanks the V, D, and J regions of immunoglobulin and T-cell receptors. Alternate splicing results in both coding and non-coding transcript variants. [provided by RefSeq, Sep 2011]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).