

## Product datasheet for **SR311962**

### **BIGM103 (SLC39A8) Human siRNA Oligo Duplex (Locus ID 64116)**

#### **Product data:**

|                            |   |
|----------------------------|---|
| <b>Product Type:</b>       | siRNA Oligo Duplexes  |
| <b>Purity:</b>             | HPLC purified   |
| <b>Quality Control:</b>    | Tested by ESI-MS  |
| <b>Sequences:</b>          | Available with shipment   |
| <b>Stability:</b>          | One year from date of shipment when stored at -20°C.  |
| <b># of transfections:</b> | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).  |
| <b>Note:</b>               | Single siRNA duplex (10nmol) can be ordered.  |
| <b>RefSeq:</b>             | <a href="#">NM_001135146</a> , <a href="#">NM_001135147</a> , <a href="#">NM_001135148</a> , <a href="#">NM_022154</a>  |
| <b>UniProt ID:</b>         | <a href="#">Q9C0K1</a>  |
| <b>Synonyms:</b>           | BIGM103; CDG2N; LZT-Hs6; PP3105; ZIP8   |
| <b>Components:</b>         | SLC39A8 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 64116)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml   |
| <b>Summary:</b>            | This gene encodes a member of the SLC39 family of solute-carrier genes, which show structural characteristics of zinc transporters. The encoded protein is glycosylated and found in the plasma membrane and mitochondria, and functions in the cellular import of zinc at the onset of inflammation. It is also thought to be the primary transporter of the toxic cation cadmium, which is found in cigarette smoke. Multiple transcript variants encoding different isoforms have been found for this gene. Additional alternatively spliced transcript variants of this gene have been described, but their full-length nature is not known. [provided by RefSeq, Oct 2008] |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).