

## Product datasheet for **SR312103**

### **RAI16 (FAM160B2) Human siRNA Oligo Duplex (Locus ID 64760)**

#### **Product data:**

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<u><a href="#">NM_022749</a></u> , <u><a href="#">NM_001354250</a></u> , <u><a href="#">NM_001354251</a></u> , <u><a href="#">NR_148727</a></u> , <u><a href="#">NR_148728</a></u> , <u><a href="#">NR_148729</a></u> , <u><a href="#">NR_148730</a></u> , <u><a href="#">NR_148731</a></u>
UniProt ID:	<u><a href="#">Q86V87</a></u>
Synonyms:	RAI16
Components:	FAM160B2 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 64760) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Able to activate MAPK/ERK and TGFB signaling pathways (PubMed:22971576). May regulate the activity of genes involved in intestinal barrier function and immunoprotective inflammation (By similarity). May play a role in cell proliferation (PubMed:22971576). [UniProtKB/Swiss-Prot Function]



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).