

Product datasheet for **SR312729**

DENND1C Human siRNA Oligo Duplex (Locus ID 79958)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_001290331 , NM_024898
UniProt ID:	Q8IV53
Synonyms:	FAM31C
Components:	DENND1C (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 79958) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	The protein encoded by this gene functions as a guanine nucleotide exchange factor for the early endosomal small GTPase RAB35, which regulates endosomal membrane trafficking and is involved in actin polymerization. The encoded protein activates RAB35 by promoting the exchange of RAB35-bound GDP for GTP. This gene may play a role in linking RAB35 activation with the clathrin machinery. [provided by RefSeq, May 2017]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).