

## Product datasheet for **SR313142**

### GSG1 Human siRNA Oligo Duplex (Locus ID 83445)

#### Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	<u><a href="#">NM_001080554</a></u> , <u><a href="#">NM_001080555</a></u> , <u><a href="#">NM_001206842</a></u> , <u><a href="#">NM_001206843</a></u> , <u><a href="#">NM_001206845</a></u> , <u><a href="#">NM_031289</a></u> , <u><a href="#">NM_153823</a></u> , <u><a href="#">NM_001367364</a></u> , <u><a href="#">NM_001367359</a></u> , <u><a href="#">NM_001367361</a></u> , <u><a href="#">NM_001367363</a></u> , <u><a href="#">NM_001367366</a></u> , <u><a href="#">NM_001368007</a></u> , <u><a href="#">NR_159953</a></u> , <u><a href="#">NM_001367358</a></u> , <u><a href="#">NM_001367360</a></u> , <u><a href="#">NM_001367362</a></u>
UniProt ID:	<u><a href="#">Q2KHT4</a></u>
Components:	GSG1 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 83445) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	May cause the redistribution of PAPOLB from the cytosol to the endoplasmic reticulum. [UniProtKB/Swiss-Prot Function]
Performance Guaranteed:	OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).



[View online »](#)