

## Product datasheet for **SR314617**

### **RAB3IP Human siRNA Oligo Duplex (Locus ID 117177)**

#### **Product data:**

|                            |  |
|----------------------------|--|
| <b>Product Type:</b>       | siRNA Oligo Duplexes   |
| <b>Purity:</b>             | HPLC purified  |
| <b>Quality Control:</b>    | Tested by ESI-MS   |
| <b>Sequences:</b>          | Available with shipment  |
| <b>Stability:</b>          | One year from date of shipment when stored at -20°C.   |
| <b># of transfections:</b> | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).   |
| <b>Note:</b>               | Single siRNA duplex (10nmol) can be ordered.   |
| <b>RefSeq:</b>             | <a href="#">NM_001024647</a> , <a href="#">NM_001278402</a> , <a href="#">NM_022456</a> , <a href="#">NM_175623</a> , <a href="#">NM_175624</a> , <a href="#">NM_175625</a> ,<br><a href="#">NM_175626</a> , <a href="#">NM_175627</a> , <a href="#">NR_103519</a> , <a href="#">NR_103520</a>   |
| <b>UniProt ID:</b>         | <a href="#">Q96QF0</a>   |
| <b>Synonyms:</b>           | RABIN3; RABIN8   |
| <b>Components:</b>         | RAB3IP (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 117177)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml  |
| <b>Summary:</b>            | Guanine nucleotide exchange factor (GEF) which may activate RAB8A and RAB8B. Promotes the exchange of GDP to GTP, converting inactive GDP-bound Rab proteins into their active GTP-bound form. Mediates the release of GDP from RAB8A and RAB8B but not from RAB3A or RAB5. Modulates actin organization and promotes polarized transport of RAB8A-specific vesicles to the cell surface. Together with RAB11A, RAB8A, the exocyst complex, PARD3, PRKCI, ANXA2, CDC42 and DNMBP promotes transcytosis of PODXL to the apical membrane initiation sites (AMIS), apical surface formation and lumenogenesis.[UniProtKB/Swiss-Prot Function] |



[View online »](#)

**Performance  
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).