

## Product datasheet for **SR314870**

### SNHG29 Human siRNA Oligo Duplex (Locus ID 125144)

#### Product data:

|                         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:           | siRNA Oligo Duplexes                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Purity:                 | HPLC purified                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Quality Control:        | Tested by ESI-MS                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Sequences:              | Available with shipment                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Stability:              | One year from date of shipment when stored at -20°C.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| # of transfections:     | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Note:                   | Single siRNA duplex (10nmol) can be ordered.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| RefSeq:                 | <a href="#">NM_152350</a> , <a href="#">NR_027158</a> , <a href="#">NR_027159</a> , <a href="#">NR_027160</a> , <a href="#">NR_027161</a> , <a href="#">NR_027162</a> , <a href="#">NR_027163</a> , <a href="#">NR_027164</a> , <a href="#">NR_027165</a> , <a href="#">NR_027166</a> , <a href="#">NR_027167</a> , <a href="#">NR_027168</a> , <a href="#">NR_027169</a> , <a href="#">NR_027170</a> , <a href="#">NR_027171</a> , <a href="#">NR_027172</a> , <a href="#">NR_027173</a> , <a href="#">NR_027174</a> , <a href="#">NR_027175</a> , <a href="#">NR_027176</a> , <a href="#">NR_027177</a> , <a href="#">NR_027178</a> , <a href="#">NR_027179</a> , <a href="#">NR_027667</a> , <a href="#">NR_045021</a> , <a href="#">NR_045022</a> , <a href="#">NR_045023</a> , <a href="#">NR_045024</a> , <a href="#">NR_045025</a> , <a href="#">NR_045026</a> , <a href="#">NR_045027</a> , <a href="#">NR_045028</a> , <a href="#">NR_045029</a> |
| Synonyms:               | FLJ25777; MGC40157; TSAP19                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Components:             | LRRC75A-AS1 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 125144)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Performance Guaranteed: | OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).



[View online »](#)