

Product datasheet for **SR314950**

DNAI3 Human siRNA Oligo Duplex (Locus ID 126820)

Product data:

Product Type:	siRNA Oligo Duplexes
Purity:	HPLC purified
Quality Control:	Tested by ESI-MS
Sequences:	Available with shipment
Stability:	One year from date of shipment when stored at -20°C.
# of transfections:	Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).
Note:	Single siRNA duplex (10nmol) can be ordered.
RefSeq:	NM_001288563 , NM_145172
UniProt ID:	Q8IWG1
Synonyms:	DIC3; NYD-SP29; WDR63
Components:	WDR63 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 126820) Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml
Summary:	Acts as a negative regulator of cell migration, invasion, and metastasis downstream of p53/TP53, through inhibition of Arp2/3 complex-mediated actin polymerization (PubMed:32128961). Via its association with the multisubunit axonemal dynein complex, is potentially involved in the regulation of cilia function (By similarity). May play a role in osteogenesis of dental tissue-derived mesenchymal stem cells (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)

**Performance
Guaranteed:**

OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at techsupport@origene.com. Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).