

## Product datasheet for **SR317062**

### FBXL22 Human siRNA Oligo Duplex (Locus ID 283807)

#### Product data:

|                            |                                                                                                                                                                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b>       | siRNA Oligo Duplexes                                                                                                                                                                                                                                    |
| <b>Purity:</b>             | HPLC purified                                                                                                                                                                                                                                           |
| <b>Quality Control:</b>    | Tested by ESI-MS                                                                                                                                                                                                                                        |
| <b>Sequences:</b>          | Available with shipment                                                                                                                                                                                                                                 |
| <b>Stability:</b>          | One year from date of shipment when stored at -20°C.                                                                                                                                                                                                    |
| <b># of transfections:</b> | Approximately 330 transfections/2nmol in 24-well plate under optimized conditions (final conc. 10 nM).                                                                                                                                                  |
| <b>Note:</b>               | Single siRNA duplex (10nmol) can be ordered.                                                                                                                                                                                                            |
| <b>RefSeq:</b>             | <a href="#">NM_203373</a> , <a href="#">NM_001367807</a> , <a href="#">NM_001367809</a> , <a href="#">NM_001367808</a>                                                                                                                                  |
| <b>UniProt ID:</b>         | <a href="#">Q6P050</a>                                                                                                                                                                                                                                  |
| <b>Synonyms:</b>           | Fbl22                                                                                                                                                                                                                                                   |
| <b>Components:</b>         | FBXL22 (Human) - 3 unique 27mer siRNA duplexes - 2 nmol each (Locus ID 283807)<br>Included - SR30004, Trilencer-27 Universal Scrambled Negative Control siRNA Duplex - 2 nmol<br>Included - SR30005, RNase free siRNA Duplex Resuspension Buffer - 2 ml |
| <b>Summary:</b>            | This gene encodes a member of the F-box protein family. This F-box protein interacts with S-phase kinase-associated protein 1A and cullin in order to form SCF complexes which function as ubiquitin ligases.[provided by RefSeq, Sep 2010]             |



**Performance Guaranteed:** OriGene guarantees that at least two of the three Dicer-Substrate duplexes in the kit will provide at least 70% or more knockdown of the target mRNA when used at 10 nM concentration by quantitative RT-PCR when the TYE-563 fluorescent transfection control duplex (cat# SR30002) indicates that >90% of the cells have been transfected and the HPRT positive control (cat# SR30003) provides 90% knockdown efficiency.

For non-conforming siRNA, requests for replacement product must be made within ninety (90) days from the date of delivery of the siRNA kit. To arrange for a free replacement with newly designed duplexes, please contact Technical Services at [techsupport@origene.com](mailto:techsupport@origene.com). Please provide your data indicating the transfection efficiency and measurement of gene expression knockdown compared to the scrambled siRNA control (quantitative RT-PCR data required).